# Silicone Finger Implant and Sizers

Instructions for Use

Mode d'emploi

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Bruksanvisning

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## Instructions for Use

### OsteoTec Silicone Finger Implant and Sizers

#### Description

The OsteoTec Silicone Finger Implant is a one-piece, flexible silicone elastomer implant which is placed in the proximal and distal intramedullary canals of either metacarpophalangeal (MCP) or proximal interphalangeal (PIP) joints after resection arthroplasty. It is designed to help restore function to hands disabled by rheumatoid, degenerative or post-traumatic arthritis.

The two intramedullary stems of the implant are joined by an integral flexible hinge which allows load distribution and causes minimal restriction to flexion/extension. The implant is not fixed into the bones, but acts as a joint spacer and becomes stabilised by the encapsulation process. The implant is available in eleven sizes to accommodate a wide range of patient anatomy and is supplied sterile.

A set of eleven implant sizers is available to determine the appropriate size of implant. They are supplied non-sterile and are not suitable for implantation.

#### General Indications

As with any joint reconstruction procedure, the following general indications should be considered before deciding whether to use this implant.

- Good patient condition
- Adequate bone stock to support the implant
- Adequate skin coverage
- Patient co-operation
- Good neurovascular status
- Potentially functional musculotendinous system
- Available post-operative care

#### Clinical Indications

This implant is intended primarily for use in resection arthroplasty of the MCP joints and, in some cases, the PIP joint, in the following cases:

In rheumatoid or post-traumatic disabilities of the MCP joint, with:

- Fixed or stiff MCP joints
- X-ray evidence of joint destruction or subluxation
- Ulnar drift, not correctable by more conservative treatment
- Contracted intrinsic and extrinsic musculature and ligament system
- Associated stiff interphalangeal joints

In rheumatoid, degenerative or post-traumatic disabilities of the PIP joint with:

- Destroyed or subluxed joint
- Stiffened joints which could not be corrected by a joint tissue release

#### Contraindications

- Infection
- Un-cooperative patient
- Inadequate condition of bone, skin or neurovascular system
- Permanently damaged tendon system
- Potentially successful conservative treatment
- Young patients with open epiphyses

#### Warnings and Precautions

- This implant should **not** be used in highly active or high-load-bearing patients.
- OsteoTec does not recommend a particular surgical technique for use with this implant. It remains the responsibility of the medical professional to evaluate and use the appropriate surgical technique based on their personal clinical training and experience.
- Only handle the implant using blunt instruments in order to avoid contamination or damage to the surface of the implant.
- Do not try to reshape the implant as this can adversely affect its structural integrity.
- Inspect the packaging before use. If the inner packaging is **not** intact, do **not** use the implant.
- Rinse the implant in sterile saline solution prior to implantation.
- The implant is for single use only – do not attempt to re-sterilise the implant as this can affect its material properties and can reduce the life of the implant.

#### Potential Complications and Adverse Side Effects

As with any joint replacement procedure, a number of complications or adverse side effects can potentially occur, these generally include but are not limited to:

- Implant loosening or fracture may necessitate revision surgery
- Excessive patient activity may cause implant wear or failure
- Infection, pain, swelling and inflammation may occur at the implant site
- Patient may have an allergic, immunological reaction, autoimmune disorder or histological response to the implant material requiring removal of the implant
- Implant movement and wear may generate wear particles that may cause or exacerbate synovitis or the formation of bone cysts in the tissues surrounding the implant.

#### Sterilisation of the Sizers

The sizers are supplied non-sterile and are not to be used as implants. The sizers must be sterilised before use according to the following method:

- Scrub using a soft-bristled brush with a solution of hot water and mild, non-oily soap to remove any surface contamination. Do not use synthetic detergents or oily soaps.
- Rinse thoroughly with distilled water.
- Sterilise using moist steam sterilisation in a high-vacuum autoclave at 134-137°C for 3.5 minutes.

Other sterilisation methods and cycles may also be suitable. However, individuals and hospitals are advised to validate whichever method they deem appropriate at their institution. Do **not** use ethylene oxide to sterilise the sizers as gas residue in the sizers can cause an adverse tissue reaction.

#### Reprocessing of single use devices

Devices labelled as single use only may not perform as intended by the manufacturer if re-used. Use of these devices cause irreversible changes to the micro and macro structure; consequently performance characteristics of the device will be sub-optimal if re-used. Re-use of a single use device may lead to:

- An increased risk of infection
- Material degradation
- Failure of the device to perform as intended
- Endotoxin reactions

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