## SkeliGraft™

**En:** Blocks and Granules

**Fr:** Bloques et granulats

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**Pl:** Bloczki i granulat

**Cz:** Bločky a granulę

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**Osteotec Limited**

9 SILVER BUSINESS PARK, AIRFIELD WAY, CHRISTCHURCH, DORSET, BH23 3TA, UK

Tel: 01202 487885 Fax: 01202 487886 **Web: www.osteotec.co.uk**

OST-ON-BG-IFU Issue 2: 050110

### INSTRUCTIONS FOR USE (En)

#### DESCRIPTION

SkeliGraft™ is a mixture having a nominal composition of 60% hydroxyapatite and 40% β-tricalcium phosphate. These materials have been subjected to extensive and successful clinical studies for many years. SkeliGraft™ is a safe and fully biocompatible material which is designed to act as an osteoconductive scaffold to support the ingrowth and fusion of adjacent viable bone when placed in an osseous environment. SkeliGraft™ will later resorb and be replaced by natural bone.

SkeliGraft™ can in many cases spare the patient the trauma of autograft harvesting and provides a sterile alternative to human or animal cadaver bone that is free of the risk of disease transmission.

#### CONDITIONS OF USE

SkeliGraft™ is for professional use only. Responsibility for proper selection of patients, for adequate training, for experience in the choice of SkeliGraft™, all aspects of the surgery, and for the choice of post-operative procedures rests entirely with the surgeon.

#### INDICATIONS

SkeliGraft™ is a synthetic osteoconductive bone substitute in block and granular form that is intended to be used as a scaffold to support the ingrowth of adjacent viable bone in defects that are not intrinsic to the stability of the bone structure. These defects may be located in the extremities, spine, pelvis or cranium.

The blocks or granules may be gently pressed into the defect by hand or shaped to more accurately fit the defect, and may also be soaked with bone marrow aspirate prior to implantation. SkeliGraft™ will act as a temporary scaffold and is **not intended to provide structural support** during the healing process. The implant is biocompatible, radio-opaque and will eventually resorb in a controlled way.

#### ADVERSE EFFECTS

A successful result may not be achieved in every case.

A secondary operation to remove or replace an implant may be necessary due to surgical error, specific medical conditions or device failure. Possible adverse effects may include but are not limited to:

- Fracture or crushing of the implant with or without generation of particulate debris due to a load being applied.
- Bone deformity or non-viable bone at the implantation site.
- Post operative complications including infection, haematoma, oedema, swelling, and fluid accumulation, tissue thinning, and other complications that are possible with surgery.
- Allergic reaction to the product.

#### PRECAUTIONS

SkeliGraft™ is for professional use only, and is only intended for use by surgeons familiar with and skilled in the techniques of bone repair and replacement.

SkeliGraft™ is not intended for load bearing applications. It is important to ensure that the area around the implantation site be secured mechanically with rigid fixation to provide structural support and maintain the implant in a static, load free environment. It is important to maximise the contact between existing bone and the implant to ensure proper bone regeneration. SkeliGraft™ must not be used to gain screw fixation.

The effect of SkeliGraft™ on patients with the following conditions is not known.

- Long term infection.
- Metabolic bone disease.
- Radiation bone therapy.
- Cardiovascular disease.

The effects of SkeliGraft™ in paediatric patients and the effects of mixing with other substances are unknown.

#### WARNINGS

Blocks and granules must be secured to prevent any potential migration of the implant and should only be used in procedures where the implant can be adequately contained.

SkeliGraft™ is opaque to x-rays and may hide areas under or above the implant on a radiograph.

#### CONTRAINDICATIONS

SkeliGraft™ is not designed or sold for any other use except as indicated. Do not use SkeliGraft™ in the presence of any contraindication.

SkeliGraft™ is contraindicated where the device is intended to provide structural support in the skeletal system and MUST NOT BE USED TO GAIN SCREW FIXATION.

Other contraindications include:

- Existing acute or chronic infections, particularly at the implantation site.
- Severe vascular or neurological disease.

- Uncontrolled diabetes.
- Severe degenerative disease.
- Hypercalcaemia, abnormal calcium metabolism.
- Inflammatory bone disease.
- Malignant tumours.
- Severely impaired renal function.

The judgement of whether an implant is suitable for a particular patient is at the surgeon’s discretion.

#### APPLICATION

**Step 1:** Open both the outer and inner pouches. Open the container (for granular product only).

**Step 2:** Implant the material. Blocks may be gently shaped to fit the defect. The implant may be first soaked with bone marrow aspirate and or blood as per normal surgical practice and may then be gently and carefully tamped into place.

**Step 3:** Secure the surgical site after implanting to prevent motion and implant migration. When excess fluid is present at the site, cauterisation, suction and application of bone wax (if needed) may be used to reduce bleeding. If the material does not fill the site correctly, remove the implant and repeat the process with a fresh amount of SkeliGraft™.

#### PACKAGING

Devices labelled sterile should only be assumed thus, so long as the manufacturer’s original packaging is intact. All SkeliGraft™ products should be received in an intact package. Damaged packages or products should not be used and should be discarded. SkeliGraft™ can NOT be re-sterilised. The expiry date is printed on the labelling. Do not use SkeliGraft™ after the expiry date.

It is recommended not to open the inner pouch until a few minutes prior to implantation.

#### STERILITY

The condition of the device is denoted by the symbol “STERILE R” Devices labelled sterile should only be assumed thus, so long as the manufacturers original packaging is intact with due regard to the expiry date shown on the product labelling.

**STERILE PRESENTATION:** Products bearing a “STERILE R” symbol have been pre-sterilised by gamma irradiation and benefit from two sterile barrier packages closest to the product. All other packaging should be regarded as for transit and storage purposes only and NOT to represent a sterile barrier.

#### RE-STERILISATION

SkeliGraft™ can NOT be re-sterilised.

#### REPROCESSING OF SINGLE USE DEVICES

Devices labelled as single use only may not perform as intended by the manufacturer if re-used. Use of these devices cause irreversible changes to the micro and macro structure; consequently performance characteristics of the device will be sub-optimal if re-used. Re-use of a single use device may lead to:

- An increased risk of infection
- Material degradation
- Failure of the device to perform as intended
- Endotoxinic reactions

#### MATERIALS

SkeliGraft™ is a mixture having a nominal composition of 60% hydroxyapatite and 40% β-tricalcium phosphate.

#### ENVIRONMENTAL STORAGE/CONDITIONS

SkeliGraft™ is intended to be stored in a dry place at room temperature. Optimal storage conditions: 15-25°C (59-77°F), 40-60% relative humidity.

#### SINGLE USE

SkeliGraft™ is intended for SINGLE USE ONLY and should under no circumstances be re-used.

#### DISPOSAL

SkeliGraft™ is environmentally friendly. No special disposal is required.

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